


Product Datasheet

p95 NBS1 (NBN) antibody (orb1318708)

Description	Mouse monoclonal antibody to p95 NBS1 (NBN)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 183-460 of human NBN (NP_002476) produced in E.coli.
Target	NBN
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	84.8 kDa
Uniprot ID	O60934
NCBI	NM_002485
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...