

Product Datasheet

DNA Polymerase iota (POLI) antibody (orb1318722)

Description	Mouse monoclonal antibody to DNA Polymerase iota (POLI)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 458-740 of human POLI (NP_009126) produced in E.coli.
Target	POLI
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	82.8 kDa
Uniprot ID	Q9UNA4
NCBI	NM_007195
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.

