


Product Datasheet

RASGRP3 antibody (orb1318740)

Description	Mouse monoclonal antibody to RASGRP3
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 380-690 of human RASGRP3 (NP_733772) produced in E.coli.
Target	RASGRP3
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	78.2 kDa
Uniprot ID	Q8IV61
NCBI	NM_170672
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...