

Product Datasheet

HSP90AB1 antibody (orb131879)

Description	Rabbit polyclonal to HSP90AB1. Heat shock protein HSP 90-beta (HSP90AB1) belongs to the heat shock p...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human ARID1A.
Preservatives	0.01% NaN ₃
Concentration	100 µg/ml BSA
Storage	-20°C
Note	For research use only
Application notes	WB: The apparent protein size on WB may be different from the calculated M.W. due to modifications. Experiment Notes: Rabbit polyclonal antibodies were produced by immunizing animals with a GST-fusion protein containing the C-terminal region of human HSP90AB1.
Clonality	Polyclonal
Purity	Affinity purification
MW	83
Uniprot ID	P08238
Dilution Range	WB:1:1000-1:3000
Expiration Date	12 months from date of receipt.