


Product Datasheet

KRAS antibody (orb1318837)

Description	Mouse monoclonal antibody to KRAS
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human KRAS (NP_203524) produced in HEK293T cell.
Target	KRAS
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	21.5 kDa
Uniprot ID	P01116
NCBI	NM_033360
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...