

Product Datasheet

RLIM antibody (orb131887)

Description	Rabbit polyclonal to RLIM. RING finger LIM domain-binding protein (RLIM) is a RING-H2 zinc finger pr...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human RLIM
Preservatives	0.01% NaN ₃
Concentration	100 µg/ml BSA
Storage	-20°C
Note	For research use only
Application notes	WB: The apparent protein size on WB may be different from the calculated M.W. due to modifications. Experiment Notes: Rabbit polyclonal antibodies were produced by immunizing animals with a GST-fusion protein containing the N-terminal region of human RLIM.
Clonality	Polyclonal
Purity	Affinity purification
MW	69
Uniprot ID	Q9NVW2
Dilution Range	WB:1:1000-1:3000
Expiration Date	12 months from date of receipt.