


Product Datasheet

HE4 (WFDC2) antibody (orb1318895)

Description	Mouse monoclonal antibody to HE4 (WFDC2)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 31-124 of human WFDC2 (NP_006094) produced in E.coli.
Target	WFDC2
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	10 kDa
Uniprot ID	Q14508
NCBI	NM_006103
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...