


Product Datasheet

ARID1A antibody pairs (orb131893)

Description	Rabbit polyclonal antibody to ARID1A.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	N-terminal or the C-terminal region of human ARID1A
Preservatives	0.01% NaN ₃
Concentration	100 µg/ml BSA
Storage	-20°C
Note	For research use only
Application notes	WB: The apparent protein size on WB may be different from the calculated M.W. due to modifications. Experiment Notes: Rabbit polyclonal antibodies were produced by immunizing animals with GST-fusion proteins containing either the N-terminal or the C-terminal region of human ARID1A.
Clonality	Polyclonal
Purity	Affinity purification
MW	242
Uniprot ID	O14497
Dilution Range	WB:1:1000-1:3000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...