


Product Datasheet

Parathyroid Hormone (PTH) antibody (orb1318966)

Description	Mouse monoclonal antibody to Parathyroid Hormone (PTH)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 32-115 of human PTH (NP_000306) produced in E.coli.
Target	PTH
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	9.4 kDa
Uniprot ID	P01270
NCBI	NM_000315
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...