

Product Datasheet

LIN37 antibody (orb131902)

Description	Rabbit polyclonal to LIN37. Protein lin-37 homolog (LIN37) is a component of the DREAM complex (also...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human LIN37
Preservatives	0.01% NaN ₃
Concentration	100 µg/ml BSA
Storage	-20°C
Note	For research use only
Application notes	WB: The apparent protein size on WB may be different from the calculated M.W. due to modifications. Experiment Notes: Rabbit polyclonal antibodies were produced by immunizing animals with a GST-fusion protein containing the N-terminal region of human LIN37.
Clonality	Polyclonal
Purity	Affinity purification
MW	39
Uniprot ID	Q96GY3
Dilution Range	WB:1:1000-1:3000
Expiration Date	12 months from date of receipt.