


Product Datasheet

MAGEA3 antibody (orb1319020)

Description	Mouse monoclonal antibody to MAGEA3
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human MAGEA3 (NP_005353) produced in HEK293T cell.
Target	MAGEA3
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	34.6 kDa
Uniprot ID	P43357
NCBI	NM_005362
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...