


Product Datasheet

CD43 (SPN) antibody (orb1319023)

Description	Mouse monoclonal antibody to CD43 (SPN)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human SPN (NP_003114) produced in HEK293T cell.
Target	SPN
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	40.1 kDa
Uniprot ID	P16150
NCBI	NM_003123
Dilution Range	WB 1:2000, IHC 1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...