


## Product Datasheet

### Ki67 (MKI67) antibody (orb1319032)

|                            |                                                                                                                          |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to Ki67 (MKI67)                                                                                |
| <b>Species/Host</b>        | Mouse                                                                                                                    |
| <b>Reactivity</b>          | Human                                                                                                                    |
| <b>Conjugation</b>         | Unconjugated                                                                                                             |
| <b>Tested Applications</b> | IHC, WB                                                                                                                  |
| <b>Immunogen</b>           | Human recombinant protein fragment corresponding to amino acids 1160-1493 of human MKI67 (NP_002408) produced in E.coli. |
| <b>Target</b>              | MKI67                                                                                                                    |
| <b>Note</b>                | For research use only                                                                                                    |
| <b>Isotype</b>             | IgG1                                                                                                                     |
| <b>Clonality</b>           | Monoclonal                                                                                                               |
| <b>Uniprot ID</b>          | <a href="#">P46013</a>                                                                                                   |
| <b>NCBI</b>                | <a href="#">NM_002417</a>                                                                                                |
| <b>Dilution Range</b>      | IHC 1:150, WB 1:200 - 1:1000                                                                                             |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                          |



Western blot analysis of extracts (10ug)...