


Product Datasheet

Hsp60 (HSPD1) antibody (orb1319044)

Description	Mouse monoclonal antibody to Hsp60 (HSPD1)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human HSPD1 (NP_955472) produced in E.coli.
Target	HSPD1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	60.9 kDa
Uniprot ID	P10809
NCBI	NM_199440
Dilution Range	WB 1:500, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...