

Product Datasheet

EGR2 antibody (orb131906)

Description	Rabbit polyclonal to EGR2. E3 SUMO-protein ligase EGR2 belongs to the EGR C2H2-type zinc-finger prot...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ChIP, IP, WB
Immunogen	C-terminal region of human EGR2
Preservatives	0.01% NaN ₃
Concentration	100 µg/ml BSA
Storage	-20°C
Note	For research use only
Application notes	WB: The apparent protein size on WB may be different from the calculated M.W. due to modifications. Experiment Notes: Rabbit polyclonal antibodies were produced by immunizing animals with a GST-fusion protein containing the C-terminal region of human EGR2.
Clonality	Polyclonal
Purity	Affinity purification
MW	50
Uniprot ID	P11161
Dilution Range	WB:1:1000-1:3000,IF/ICC:1:100-1:200
Expiration Date	12 months from date of receipt.