

Product Datasheet

ASCL3 antibody (orb131909)

Description	Rabbit polyclonal to ASCL3. Achaete-scute homolog 3 (ASCL3) is a class B basic helix-loop-helix (bHL...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	C-terminal region of mouse ASCL3
Preservatives	0.01% NaN ₃
Concentration	100 µg/ml BSA
Storage	-20°C
Note	For research use only
Application notes	WB: The apparent protein size on WB may be different from the calculated M.W. due to modifications. Experiment Notes: Rabbit polyclonal antibodies were produced by immunizing animals with a GST-fusion protein containing the C-terminal region of mouse ASCL3.
Clonality	Polyclonal
Purity	Affinity purification
MW	21
Uniprot ID	Q9JJR7
Dilution Range	WB:1:1000-1:3000,IF/ICC:IP:1:100-1:200
Expiration Date	12 months from date of receipt.