


Product Datasheet

ING2 antibody (orb1319105)

Description	Mouse monoclonal antibody to ING2
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human ING2 (NP_001555) produced in E.coli.
Target	ING2
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	32.6 kDa
Uniprot ID	Q9H160
NCBI	NM_001564
Dilution Range	WB 1:500~2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...