


Product Datasheet

Ki67 (MKI67) antibody (orb1319114)

Description	Mouse monoclonal antibody to Ki67 (MKI67)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 2924-3256 of human MKI67 (NP_002408) produced in E.coli.
Target	MKI67
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P46013
NCBI	NM_002417
Dilution Range	IHC 1:50~500, WB 1:200 - 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts (10ug)...