


Product Datasheet

TIA1 antibody (orb1319117)

Description	Mouse monoclonal antibody to TIA1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human TIA1 (NP_071505) produced in HEK293T cell.
Target	TIA1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P31483
NCBI	NM_022173
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...