


Product Datasheet

c-Jun (JUN) antibody (orb1319138)

Description	Mouse monoclonal antibody to c-Jun (JUN)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human JUN (NP_002219) produced in E.coli.
Target	JUN
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	35.5 kDa
Uniprot ID	P05412
NCBI	NM_002228
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...