


Product Datasheet

CCDC93 antibody (orb1319170)

Description	Mouse monoclonal antibody to CCDC93
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 293-631 of human CCDC93 (NP_061917) produced in E.coli.
Target	CCDC93
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	73 kDa
Uniprot ID	Q567U6
NCBI	NM_019044
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...