

Product Datasheet

PRR11 antibody (orb1319183)

Description	Mouse monoclonal antibody to PRR11
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 125-360 of human PRR11 (NP_060774) produced in E.coli.
Target	PRR11
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	39.9 kDa
Uniprot ID	Q96HE9
NCBI	NM_018304
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.

