


Product Datasheet

ZFAND5 antibody (orb1319220)

Description	Mouse monoclonal antibody to ZFAND5
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human ZFAND5 (NP_005998) produced in E.coli.
Target	ZFAND5
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	23 kDa
Uniprot ID	O76080
NCBI	NM_006007
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...