


Product Datasheet

MAGEB1 antibody (orb1319229)

Description	Mouse monoclonal antibody to MAGEB1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human MAGEB1 (NP_796379) produced in E.coli.
Target	MAGEB1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	38.9 kDa
Uniprot ID	P43366
NCBI	NM_177404
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...