


Product Datasheet

A4GNT antibody (orb1319281)

Description	Mouse monoclonal antibody to A4GNT
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 121-340 of human A4GNT (NP_057245) produced in E.coli.
Target	A4GNT
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	39.3 kDa
Uniprot ID	Q9UNA3
NCBI	NM_016161
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...