


Product Datasheet

APE1 (APEX1) antibody (orb1319319)

Description	Mouse monoclonal antibody to APE1 (APEX1)
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-242 of human APEX1 (NP_001632) produced in E.coli.
Target	APEX1
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	35.4 kDa
Uniprot ID	P27695
NCBI	NM_001641
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...