


Product Datasheet

APE1 (APEX1) antibody (orb1319320)

Description	Mouse monoclonal antibody to APE1 (APEX1)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-242 of human APEX1 (NP_001632) produced in E.coli.
Target	APEX1
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	35.4 kDa
Uniprot ID	P27695
NCBI	NM_001641
Dilution Range	WB 1:1000, IHC 1:150, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Flow
cytometric
Flow
cytometric
HEK293T ...