

Product Datasheet

TLE 1 (TLE1) antibody (orb1319340)

Description	Mouse monoclonal antibody to TLE 1 (TLE1)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 180-460 of human TLE1 (NP_005068) produced in E.coli.
Target	TLE1
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	83 kDa
Uniprot ID	Q04724
NCBI	NM_005077
Dilution Range	WB 1:2000, IHC 1:150, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.

