


Product Datasheet

FBXO42 antibody (orb1319347)

Description	Mouse monoclonal antibody to FBXO42
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 356-717 of human FBXO42 (NP_061867) produced in E.coli.
Target	FBXO42
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	77.7 kDa
Uniprot ID	Q6P3S6
NCBI	NM_018994
Dilution Range	WB 1:2000, IHC 1:150, FLOW 1:100
Expiration Date	12 months from date of receipt.



Flow
cytometric
Flow
cytometric
HEK293T ...