


Product Datasheet

THAP6 antibody (orb1319366)

Description	Mouse monoclonal antibody to THAP6
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human THAP6 (NP_653322) produced in E.coli.
Target	THAP6
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	25.5 kDa
Uniprot ID	Q8TBBO
NCBI	NM_144721
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...