


Product Datasheet

MRPS27 antibody (orb1319389)

Description	Mouse monoclonal antibody to MRPS27
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 151-414 of human MRPS27 (NP_055899) produced in E.coli.
Target	MRPS27
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	47.4 kDa
Uniprot ID	Q92552
NCBI	NM_015084
Dilution Range	WB 1:200, IHC 1:150, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...