


Product Datasheet

SLUG (SNAI2) antibody (orb1319409)

Description	Mouse monoclonal antibody to SLUG (SNAI2)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human SNAI2 (NP_003059) produced in E.coli.
Target	SNAI2
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	29.8 kDa
Uniprot ID	O43623
NCBI	NM_003068
Dilution Range	WB 1:2000, IHC 1:150, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...