


Product Datasheet

TRIP (TRAIP) antibody (orb1319507)

Description	Mouse monoclonal antibody to TRIP (TRAIP)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 67-313 of human TRAIP (NP_005870) produced in E.coli.
Target	TRAIP
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	53.1 kDa
Uniprot ID	Q9BWF2
NCBI	NM_005879
Dilution Range	WB 1:2000, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...