


Product Datasheet

TRIM9 antibody (orb1319536)

Description	Mouse monoclonal antibody to TRIM9
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 284-669 of human TRIM9 (NP_443210) produced in E.coli.
Target	TRIM9
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	61.2 kDa
Uniprot ID	Q9C026
NCBI	NM_052978
Dilution Range	WB 1:2000, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...