


Product Datasheet

cIAP2 (BIRC3) antibody (orb1319549)

Description	Mouse monoclonal antibody to cIAP2 (BIRC3)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 255-529 of human BIRC3 (NP_001156) produced in E.coli.
Target	BIRC3
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	68.2 kDa
Uniprot ID	Q13489
NCBI	NM_001165
Dilution Range	WB 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...