


Product Datasheet

DPPA4 antibody (orb1319552)

Description	Mouse monoclonal antibody to DPPA4
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 9-254 of human DPPA4 (NP_060659) produced in E.coli.
Target	DPPA4
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	33.4 kDa
Uniprot ID	Q7L190
NCBI	NM_018189
Dilution Range	WB 1:500~1000, IF 1:100
Expiration Date	12 months from date of receipt.



immunofluorescent
staining of COS7
cells...