

Product Datasheet

SARS-CoV-2 N Protein antibody (orb1319602)

Description	Rabbit monoclonal antibody to SARS-CoV-2 N Protein
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Full length recombinant SARS-CoV-2 N protein produced in HEK293.
Target	N Protein
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	47kDa
NCBI	YP_009724397
Dilution Range	ELISA 1:5000-10000
Expiration Date	12 months from date of receipt.