

Product Datasheet

S Protein antibody (orb1319618)

Description	Rabbit monoclonal antibody to S Protein
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant fusion protein of SARS-CoV-2 Spike RBD Protein and Virion Surface Domain of Membrane Protein (M-protein) produced in HEK293
Target	S Protein
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	34kDa
NCBI	YP_009724390
Dilution Range	ELISA 1:5000-10000
Expiration Date	12 months from date of receipt.