

## Product Datasheet

### SARS-CoV-2 N Protein antibody (orb1319625)

|                            |                                                                  |
|----------------------------|------------------------------------------------------------------|
| <b>Description</b>         | Rabbit monoclonal antibody to SARS-CoV-2 N Protein               |
| <b>Species/Host</b>        | Rabbit                                                           |
| <b>Reactivity</b>          | Virus                                                            |
| <b>Conjugation</b>         | Unconjugated                                                     |
| <b>Tested Applications</b> | ELISA                                                            |
| <b>Immunogen</b>           | Full length recombinant SARS-Cov-2 N protein produced in E.coli. |
| <b>Target</b>              | N Protein                                                        |
| <b>Note</b>                | For research use only                                            |
| <b>Isotype</b>             | IgG                                                              |
| <b>Clonality</b>           | Monoclonal                                                       |
| <b>MW</b>                  | 47kDa                                                            |
| <b>NCBI</b>                | <b>YP_009724397</b>                                              |
| <b>Dilution Range</b>      | ELISA 1:5000-10000                                               |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                  |