


Product Datasheet

CD20 (MS4A1) antibody (orb1319705)

Description	Mouse monoclonal antibody to CD20 (MS4A1)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	A full length cDNA gene delivered by DNA immunization
Target	MS4A1
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	P11836
NCBI	NM_021950
Dilution Range	Use at 0.1-1ug/ml. Flow , CDC
Expiration Date	12 months from date of receipt.



Flow cytometry on CD20 transfected HEK28...