

Product Datasheet

ASS1 antibody (orb1319724)

Description	Rabbit monoclonal antibody to ASS1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	MS
Immunogen	Synthetic peptide (the amino acid sequence is considered to be commercially sensitive) within Human ASS1 (NP_446464). The exact sequence is proprietary.
Target	ASS1
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	46.3 kDa
Uniprot ID	P00966
NCBI	NM_054012
Expiration Date	12 months from date of receipt.