


Product Datasheet

Dendra2 antibody (orb1319748)

Description	Rabbit monoclonal antibody to Dendra2
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	This antibody is raised against a recombinant protein of Dendra2 expressed in HEK293 cells
Target	Dendra2
Preservatives	Tissue culture buffer containing 1% rAlbumin and 0.02% Sodium azide
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	27.5 kDa
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of Dendra2 purified...