


Product Datasheet

PSCD4 (CYTH4) antibody (orb1319760)

Description	Rabbit polyclonal antibody to PSCD4 (CYTH4)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	DNA immunization. This antibody is specific for the N Terminus Region of the target protein.
Target	CYTH4
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UIA0
NCBI	NM_013385
Dilution Range	WB 1:20000,ELISA 1:100-1:2000
Expiration Date	12 months from date of receipt.



ELISA using
PSCD4
(CYTH4)
antibody