


Product Datasheet

PHLPP2 antibody (orb1319761)

Description	Rabbit polyclonal antibody to PHLPP2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	DNA immunization. This antibody is specific for the C Terminus Region of the target protein.
Target	PHLPP2
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6ZVD8
NCBI	NM_015020
Dilution Range	WB 1:20000, ELISA 1:100-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...