


Product Datasheet

KMT2A antibody (orb1319768)

Description	Rabbit polyclonal antibody to KMT2A
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	DNA immunization. This antibody is specific for the Middle Region of the target protein.
Target	KMT2A
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q03164
NCBI	NM_001197104
Dilution Range	ELISA: 1:100-1:2000; IHC: 1:150
Expiration Date	12 months from date of receipt.



ELISA using
KMT2A
antibody