

Product Datasheet

EV11 (MECOM) antibody (orb1319771)

Description	Rabbit polyclonal antibody to EV11 (MECOM)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	DNA immunization. This antibody is specific for the N Terminus Region of the target protein.
Target	MECOM
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q03112
NCBI	NM_001105077
Dilution Range	WB: 1:5000-1:20000; ELISA: 1:100-1:2000; IHC: 1:10-1:2000; IHC-P 1:250-1:2000
Expiration Date	12 months from date of receipt.



ELISA using
EV11
(MECOM)
antibody