


Product Datasheet

TBC1D8 antibody (orb1319815)

Description	Rabbit polyclonal antibody to TBC1D8
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC
Immunogen	DNA immunization. This antibody is specific for the N Terminus Region of the target protein.
Target	TBC1D8
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O95759
NCBI	NM_001102426
Dilution Range	ELISA: 1:100-1:2000; IHC: 1:10-1:2000; IHC-P 1:10-1:2000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human ga...