


Product Datasheet

CLNS1A antibody (orb1319817)

Description	Rabbit polyclonal antibody to CLNS1A
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	DNA immunization. This antibody is specific for the N Terminus Region of the target protein.
Target	CLNS1A
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P54105
NCBI	NM_001293
Dilution Range	WB 1:5000~20000, ELISA 1:100-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...