

Product Datasheet

Osteoprotegerin (TNFRSF11B) antibody (orb1319844)

Description Rabbit polyclonal antibody to Osteoprotegerin (TNFRSF11B)

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications ELISA

Immunogen DNA immunization. This antibody is specific for the C Terminus Region of the target protein.

Target TNFRSF11B

Note For research use only

Isotype IgG


Clonality Polyclonal

Uniprot ID [O00300](#)

NCBI [NM_002546](#)

Dilution Range WB: 1:5000-1:20000; ELISA: 1:100-1:2000

Expiration Date 12 months from date of receipt.



ELISA using
Osteoprotegerin
(TNFRSF11B) ...