


Product Datasheet

PCK1 antibody (orb1319845)

Description	Rabbit polyclonal antibody to PCK1
Species/Host	Rabbit
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	DNA immunization. This antibody is specific for the Middle Region of the target protein.
Target	PCK1
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P35558
NCBI	NM_002591
Dilution Range	WB 1:5000~20000, IHC 1:150, ELISA 1:100-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...